

# ENDOCEPT™

**Multi-function enzymatic detergent for manual and automated reprocessing of surgical instruments and flexible endoscopes.**



**Superior cleaning performance**

Ensuring reliable results every time.



**Premium formulation**

Removes organic proteins and carbohydrates delivering exceptional results.



**Highly concentrated formula**

Means more uses per bottle – less waste, more value.



**Passes a range of commercially available soil tests**

To ensure cleaning across a wide range of applications.



**Phosphate-free formula**

Ensures powerful cleaning without compromising sustainability.

Endocept™ is a multi-function enzymatic detergent for manual and automated reprocessing of surgical instruments and flexible endoscopes.

Endocept™ is a low foaming, pH neutral, non-corrosive and free-rinsing formula compatible with reusable surgical instruments and flexible endoscopes. Endocept™ is effective at penetrating and breaking down the most common soils found on medical devices including protein, fat and carbohydrate-based soils.

With excellent materials compatibility, Endocept™ is an advanced multi-enzyme system designed for safe and effective use across a wide range of medical instruments and devices. Its formulation is engineered to be non-corrosive and gentle on delicate surfaces including stainless steel, plastics, and complex device components, helping to preserve the integrity and longevity of reusable instruments. Endocept™ delivers proven performance producing consistent results in validated cleaning tests.

Endocept™ rapidly removes protein-rich biological soils and bioburden and passes most commercial soil tests.

## Designed for universal applications

- Manual cleaning
- Ultrasonic cleaning
- Automatic Endoscope Reprocessors

## Advanced Technology

- Utilizes the latest enzyme technology specific for surgical soils
- Optimized surfactants for removal of proteins
- Effective in both hard and soft water across a range of temperatures

For more information contact Whiteley

Toll Free Phone: (888) 655-9194  
Email: [contact@whiteleyusa.com](mailto:contact@whiteleyusa.com)  
[www.whiteleyusa.com](http://www.whiteleyusa.com)



Endocept™ is specifically formulated with a blend of high-quality medical-grade enzymes to effectively remove a wide range of debris and soils and provides excellent cleaning efficacy.

The first and often vital step in any process of sterilisation is the removal of organic soils. Meticulous cleaning prior to disinfecting or sterilisation is essential. If blood or other organic matter is not entirely removed prior to disinfection or sterilisation, sterility of the equipment may be compromised. Failure to remove organic soiling prior to disinfection or sterilisation may result in the transmission of infectious organisms. Endocept™ acts quickly by penetrating and degrading dried blood and other proteinaceous matter that is difficult to remove by traditional cleaning methods.

## Benefits

- Low-foaming, crystal clear in-use solution for manual cleaning allows visibility of instruments in the sink
- Neutral pH, non-corrosive and free-rinsing formulation is compatible with surgical instruments and flexible endoscopes
- Effectively removes hard-to-see soil residues that can remain after cleaning
- Versatile and multi-purpose – can be used as a pre-soak, manual cleaner, ultrasonic detergent, and flexible endoscope detergent in hard & soft water
- Concentrated formula so less product is required to maintain performance
- Phosphate-free formula ensures powerful cleaning without compromising sustainability

For more information contact Whiteley

Toll Free Phone: (888) 655-9194  
Email: [contact@whiteleyusa.com](mailto:contact@whiteleyusa.com)  
[www.whiteleyusa.com](http://www.whiteleyusa.com)